



Clinical Biochemistry

Lab. 5

Measurement of Serum Creatinine

Prepared by :

Darya Shorsh Hamad

Mcs. in Clinical Biochemistry

E-mail: darya.shorsh@sulicihan.edu.krd

Introduction



Clinical Significance

- ❑ **Creatinine** is synthesized in the body at a **fairly constant rate** from **creatine**.
- ❑ Creatinine in the blood is then removed by filtration through the **glomeruli of the kidney** for excretion in the urine.
- ❑ the *creatinine clearance (CC) test* is one of the **most sensitive tests** to **diagnose renal function** especially the **glomerular filtration rate (GFR)**.
- ❑ **Elevated levels of creatinine** in serum are usually associated with **renal diseases**, especially those related to **GFR such as glomerular nephritis**.


Quantitative measurement of Creatinine



Aim: Quantitative measurement of Creatinine in the serum

Principle

This procedure is based upon a modification of the original picrate reaction (Jaffe). Creatinine under alkaline conditions reacts with picrate ions forming a **reddish complex**. The formation rate of the complex measured through the increase of absorbance in a prefixed interval of time is proportional to the concentration of creatinine in the sample.



Reagent Composition



R1: Picric acid

R2: Alkaline Reagent: NaOH

R3: Standard: 2mg/dL

REAGENT PREPARATION

Working reagent. Mix 1 volume of **R1** + 1 volume of **R2**. Stable for 1 week at room temperature, stored **tightly closed and protected from light**.

Samples



- ✓ Serum or heparinized plasma, and urine.
- ✓ Creatinine in serum or plasma is stable up to 24 hours at 2-8°C.
Freeze for longer storage.
- ✓ Creatinine from random samples of urine is stable for 4 days at 2-8°C. Freeze for longer storage.
- ✓ The 24-hour urine samples for the Clearance Test should be collected on a **preservative (fluoride thymol)** and **immediately refrigerated.**

Procedure



1. Adjust the instrument to zero with distilled water

	Blank	Standard	Sample
Working Reagent	1ml	1ml	1ml
Standard Solution		100 μ l	
Distilled Water	100 μ l		
Serum			100 μ l

1. Mix and start stop watch
2. Read the absorbance for each tubes at 592 nm **(A1) after 30 seconds** and **after 90 seconds (A2)** of the sample addition
3. Calculate ΔA **(A2-A1)**

Calculation



$$\text{Creatinine (mg/dL)} = \frac{(\Delta A) \text{ Sample} - (\Delta A) \text{ Blank}}{(\Delta A) \text{ Standard} - (\Delta A) \text{ Blank}} \times \text{Standard Concentration}$$

Standard Concentration: **2mg/dl**

Reference Intervals:

Male: 0.7 - 1.4 mg/dL

Female: 0.6 - 1.1mg/dL

Calculation of Creatinine Clearance



$$\text{mL/min} = \frac{\text{mg creatinine/ dL URINE} \times \text{mL 24-h}}{\text{mg creatinine/ dL SERUM} \times 1440 \text{ min}}$$

Cockcroft-Gault Formula for Estimating Creatinine Clearance

$$\text{CrCl (mL/min)} = \frac{(140 - \text{age}) \times \text{Lean Body Weight (kg)}}{\text{Serum Creatinine (mg/dL)} \times 72} \quad (\times 0.85 \text{ if female})$$

Reference Intervals:

Clearance Test

Men	97 - 137 mL/min
Women	88 - 128 mL/min