



Clinical Biochemistry

Lab. 10

Measurement of Triglyceride (TG)

Prepared by :

Darya Shorsh Hamad

Mcs. in Clinical Biochemistry



E-mail: darya.shorsh@sulicihan.edu.krd

Introduction



Blood lipids

Concept: All the lipids contained in plasma, including fat, phospholipids, cholesterol, cholesterol ester and fatty acid.



Blood lipid exist and transport in the form of **lipoprotein**.

Clinical Significance of TG Measurement



The measurement of TG concentration in the blood is important to **diagnosis and following up of Hyperlipidemia**

❖ Cause of Hyperlipidemia :

- Primary: Genetic Origin
- Secondary:
 - Diabetes Mellitus
 - Hypothyroidism

❖ Effects of Hyperlipidemia :

- Hypertriglyceridemia increases the risk for **pancreatitis**.
- Hypertriglyceridemia is associated with the following clinical findings: **xanthoma**, **Lipemia retinalis (LR)**, **hepatomegaly**, **splenomegaly** and **decreased HDL-cholesterol**.

Friedewald equation

$$\text{VLDL-c} = \text{Triglycerides}/5$$


$$\text{LDL-c} = \text{Total cholesterol} - (\text{HDL-c} + \text{VLDL-c})$$

Procedure



PRINCIPLE (4) (5)

Fossati and Prencipe method associated with Trinder reaction.
Reaction scheme is as follows:



The absorbance of the coloured complex (quinoneimine), proportional to the amount of triglycerides in the specimen, is measured at 500 nm.

Reagents



REAGENTS

Vial R1

BUFFER

PIPES	100	mmol/L
Magnesium chloride	9.8	mmol/L
Chloro-4-phenol	3.5	mmol/L
Preservative		

Vial R2

ENZYMES

Lipase	≥ 1000	IU/L
Peroxydase (POD)	≥ 1700	IU/L
Glycerol 3 phosphate oxydase (GPO)	≥ 3000	IU/L
Glycerol Kinase (GK)	≥ 660	IU/L
4 - Amino – antipyrine (PAP)	0.5	mmol/L
Adenosine triphosphate Na (ATP)	1.3	mmol/L

Vial R3

STANDARD

Glycerol	2.28	mmol/L
Equivalent to trioleine or triglycerides	200 mg/dL	(2.28 mmol/L)

Procedure



REAGENTS PREPARATION

Vial R2: Use a non-sharp instrument to remove aluminium cap.

Add promptly the contents of vial R2 (Enzymes), into vial R1 (Buffer).

Mix gently and wait for complete dissolution before using reagent (approximately 2 minutes).

Procedure



MANUAL PROCEDURE

Let stand reagent and specimens at room temperature.

Pipette into well identified test tubes:	Blank	Standard	Assay
Reagent	1 mL	1 mL	1 mL
Demineralised water	10 μ L		
Standard		10 μ L	
Specimen			10 μ L

Mix. Let stand for 5 minutes at 37°C or 10 minutes at room temperature.
Record absorbance at 500 nm (480-520) against reagent blank.
Reaction is stable for 1 hour.

Calculation



CALCULATION

Calculate the result as follows:

$$\text{Result} = \frac{\text{Abs (Assay)}}{\text{Abs (Standard)}} \times \text{Standard concentration}$$

Ranges for lipid profile Test



LIPID PROFILE

	DESIRABLE	BORDERLINE	HIGH RISK
Cholesterol	<200 mg/dl	200-239 mg/dl	240 mg/dl
Triglycerides	<150 mg/dl	150-199 mg/dl	200-499 mg/dl
HDL cholesterol	60 mg/dl	35-45 mg/dl	<35 mg/dl
LDL cholesterol	60-130 mg/dl	130-159 mg/dl	160-189 mg/dl
Cholesterol/ HDL ratio	4.0	5.0	6.0