



# Clinical Biochemistry

## Lab 3

## Measurement of Blood Glucose

Prepared by :

**Darya Shorsh Hamad**

Mcs. in Clinical Biochemistry

E-mail: [darya.shorsh@sulicihan.edu.krd](mailto:darya.shorsh@sulicihan.edu.krd)

# Introduction



- ❑ **Glucose (C<sub>6</sub>H<sub>12</sub>O<sub>6</sub>)** is a ubiquitous **fuel molecule in biology**.
- ❑ It is oxidized through a series of enzyme-catalyzed reactions to **form carbon dioxide and water**, yielding the universal energy molecule **ATP**.
- ❑ Due to its importance in metabolism, glucose level is a key diagnostic parameter for many **metabolic disorders**.
- ❑ Increased glucose levels of glucose (**Hyperglycemia**) have been associated with **diabetes mellitus, hyperactivity of thyroid, pituitary and adrenal glands**.
- ❑ Decreased levels of glucose (**Hypoglycemia**) are found in **insulin secreting tumors, myxedema, hypopituitarism and hypoadrenalism**.

## Measuring blood glucose is used for the diagnosis of:



- Carbohydrate metabolism disorders
- Monitoring of treatment in diabetes mellitus
- Neonatal hypoglycemia
- Idiopathic hypoglycemia
- Insulinoma.

# Glucose Measurements



## Glucose Measurement Methods

- ❑ Glucose is measured in **whole blood, plasma, serum, cerebrospinal fluid, pleural fluid, and urine** for a variety of diagnostic and management purposes.
- ❑ The standard clinical laboratory analysis of glucose is performed on **plasma or serum** derived from a **phlebotomy specimen**.
- ❑ Most measurements of glucose employ enzymatic methods:
  1. **Glucose Oxidase**
  2. **Hexokinase**
  3. **rarely, Glucose Dehydrogenase**

# Specimen Considerations



- ✓ Glycolysis causes plasma glucose to decline over time while the plasma or serum is in contact with cells.
- ✓ The decline can be 10 mg/dL (average 5%–7%) per hour.
- ✓ A specimen is appropriate for glucose analysis if serum or plasma is separated from the cells **within 30 minutes.**
- ✓ Until separation has been performed, it is recommended that specimens be kept on **ice.**
- ✓ Whole blood tends to give approximately **10% to 15%** lower glucose readings than plasma.
- ✓ Fluoride and iodoacetate have been used as inhibitors of glycolysis to preserve blood that cannot be separated rapidly.
- ✓ **Note: Whole blood and hemolysis are not recommended for use as a sample.**

# Blood Glucose Tests




- ❑ **Random Blood Test (RBS):** measures blood glucose regardless of when you last ate.
- ❑ **Fasting Blood Sugar (FBS):** measures blood glucose after you have not eaten for at least 8 hours. It often is the **first test done to check and diagnosis of diabetes.**
- ❑ **Glucose tolerance test (GTT):** This measures your body's ability to use sugar after drinking a standard amount in a sugary drink (75 gr). The glucose tolerance test can **be used to screen for type 2 diabetes.** evaluating an individual's apparent insulin sensitivity and insulin resistance.
- ❑ **Two Hours Postprandial Plasma glucose (2HPP):**  
Measures blood glucose exactly 2 hours after you eat a meal.it is used for Screening and monitoring of diabetes.
- ❑ **Glycosylated Hemoglobin (HA1C) :** This measures your average blood sugar level over the last 2 to 3 months

# Glycated Hemoglobin (HbA1C)



- **Glycated hemoglobin**, also known as HbA1c, glycohemoglobin, hemoglobin A1c, A1C, is a form of hemoglobin (Hb) that is chemically linked to glucose.
- It is formed in a non-enzymatic glycation pathway by **hemoglobin's exposure to plasma glucose**.
- **Normal levels of glucose produce a normal amount of glycated hemoglobin.**



- ✓ As the average amount of plasma glucose increases, the fraction of glycated hemoglobin increases in a predictable way.

# HbA1C Test



A hemoglobin A1C (HbA1C) test is a blood test that shows what your average blood sugar (glucose) level was over the past **two to three months**.

## Clinical Importance of HbA1C test

- 1. Diagnose prediabetes.** If you have prediabetes, you have a higher risk of developing diabetes and cardiovascular disease.
- 2. Diagnose type 1 and type 2 diabetes.** Along with other test like FBS is used to diagnosis of Diabetes.
- 3. Monitor your diabetes treatment plan.** The result of an initial A1C test also helps establish your baseline A1C level. The test is then repeated regularly to monitor your diabetes treatment plan.

# Diagnosis of Diabetes Mellitus



To diagnose diabetes, **two abnormal results** from following tests are required:

- (1) Fasting Plasma Glucose or FBS test:** level of **126 mg/dL (7.0 mmol/L)** or **higher** (the glucose level should be obtained after an 8-hour fast),
- (2) Oral Glucose Tolerance Test (OGTT) test :** a **plasma glucose level  $\geq 200$  mg/dL (11.1 mmol/L)** **2 hours after a glucose load** (75-g oral glucose tolerance test).
- (3) Random Blood Sugar (RBS) test :** a **Random Glucose  $\geq 200$  mg/dL (11.1 mmol/L)** with **symptoms of hyperglycemia**
- (4) Hemoglobin A1c (HbA1c) test:  $\geq 6.5\%$ .**

Tests can be performed on two different days or 2 different tests can be performed on the same day.

# Diagnosis of Diabetes Mellitus



|                            | FASTING PLASMA GLUCOSE |         | 2-HOUR PLASMA GLUCOSE LEVEL (AFTER 75 G GLUCOSE LOAD) |          | HBA <sub>1C</sub> |
|----------------------------|------------------------|---------|---|----------|-------------------|
|                            | mg/dL                  | mmol/L  | mg/dL   | mmol/L   | %                 |
| Normal                     | <100                   | <5.6    | <140  | <7.8     |                   |
| <b>Prediabetes</b>         |                        |         |   |          | <b>5.7–6.4</b>    |
| Impaired fasting glucose   | 100–125                | 5.6–6.9 |   |          |                   |
| Impaired glucose tolerance |                        |         | 140–199   | 7.8–11.0 |                   |
| Diabetes mellitus          | ≥126                   | ≥7.0    | ≥200  | ≥11.1    | ≥6.5              |

**Prediabetes** designates conditions in which **glucose homeostasis is abnormal**, but serum glucose levels are **not high enough to be classified as diabetes**. This group includes individuals with **impaired fasting glucose and impaired glucose tolerance**.

✓ They are also at increased risk for cardiovascular and cerebrovascular diseases.

# Diagnosis of Gestational Diabetes Mellitus



Perform a 75-g OGTT, with plasma glucose measurement fasting and at 1 and 2 h, at 24–28 weeks of gestation in women not previously diagnosed with overt diabetes.

The OGTT should be performed in **the morning after an overnight fast of at least 8 h**.

The diagnosis of GDM is made when any of the following plasma glucose values are exceeded:


- **Fasting:  $\geq 92$  mg/dL (5.1 mmol/L)**
- **1 h:  $\geq 180$  mg/dL (10.0 mmol/L)**
- **2 h:  $\geq 153$  mg/dL (8.5 mmol/L)**

# Methods of Glucose Analysis



- **Method 1: Benedict's:** (qualitative, semiquantitative):
- **Method 2: Alkaline ferricyanide** (quantitative): involves the reduction of yellow ferricyanide,  $\text{Fe}(\text{CN})_6$ , to colorless ferrocyanide,  $\text{Fe}(\text{CN})_6$ , by glucose in alkaline conditions.
- **Method 3: o-Toluidine** (quantitative): The o-toluidine reaction is based on the ability of many aromatic amines in acid solutions to condense with the aldehyde group of glucose to form glucosamines.  $\text{O-toluidine} + \text{glucose (aldehyde)} \xrightarrow{\text{heat \& acidity}} \text{glucosamine (colored)}$
- **Enzymatic methods: Glucose oxidase coupled reaction** (quantitative): One of the most frequently used specific glucose methods uses two coupled enzyme reactions

# Glucose oxidase coupled reaction



The intensity of the color formed is proportional to the glucose concentration in the sample. It is determined by measuring the increase in absorbance at **500-550 nm**.



# Introduction




Set the slide (Dry Slide Reagents)



Set the sample



Press START




# Fujifilm DRI-CHEM SLIDE methods



## 1. Colorimetric method slide:

(Enzymes, General chemistry, and Immunology)

This multilayered slide is composed of dry chemical ingredients needed for the reaction and other functional materials. It quantifies enzymes and chemicals using colorimetric method.




# Fujifilm DRI-CHEM SLIDE



## 2. Potentiometric method slide (Electrolytes)

Each slide comes with an ion selective film electrode for each of Na, K, and Cl.

Slides quantify electrolytes in the sample by a potentiometric method.



# Digital Glucometer



## Glucometer (or glucose meter)

Glucometer is a medical tool that is used for measuring the approximate level of glucose in the blood.

## Random Blood Glucose Test:

**Normal: 70-140 mg/ dl**

**Impaired: 140 – 200 mg/dl**

**Diabetes: > 200 mg/dl**



### OPERATION STEPS



# Working principle of a Glucometer:

---

Most of the glucometers are based on **electrochemical technology**.

## Glucometer test strips:

In each glucometer test strip which are to be used for determining glucose in the blood, contains an **enzyme called glucose oxidase**. This **enzyme then reacts** with the glucose in the blood sample and creates an acid called **gluconic acid**.

The gluconic acid thus formed then reacts with another chemical in the testing strip called **ferricyanide**.

The ferricyanide and the gluconic acid then combines with each other and forms **ferrocyanide**.

iii) As soon as the ferrocyanide has been formed the device (i.e., glucometer) **runs an electronic current through the blood sample on the strip**.

**This current thus generated is able to read the ferrocyanide and identify the amount of glucose present in the blood sample on the testing strip.**

That number is then **displayed on the screen of the glucometer**.