

Practical Microbiology

**Cihan University
Medical Laboratory Analysis**

Lab 9: Microbial counting methods (viable, direct & spectrophotometric)

Lec: Karokh Ali Khdir

MSc: Microbiology

Second stage (1st semester)

2023-2024


Objective

Each student should be able to:

- ✓ Count the bacterial cells in different samples.
- ✓ Know the methods in bacterial counting:
 - ✓ Direct microscopic count method.
 - ✓ Standard plate count method.
 - ✓ Spectrophotometric method.

Bacterial counting


- Is the enumeration of bacteria especially in the Environmental samples, clinical samples, dairy microbiology, food microbiology, and water microbiology.
- Bacterial growth may be determined based on **increase in the number of cells**.
- Enumeration methods may yield either **total counts or viable counts**.
- **Total count** is a count of cells including **dead and live** cells.
- **Viable count** is a count of only those cells that are alive in the sample.



Bacterial counting

➤ The time it takes for a population of bacterial cells to double in number is called the **generation time**.

Generation time differs significantly from one bacterium to another.



Bacterial counting

➤ Different samples could be used in bacterial counting:

- ✓ Clinical samples or body fluids (urine, CSF).
- ✓ Broth culture of bacteria.
- ✓ Environmental samples (soil).
- ✓ Water (sewage, tap water, drinking water).
- ✓ Foods (vegetables, canned food).
- ✓ Drinkable items (fruit juice, soft drink, milk, etc...).



Methods in bacterial counting

Direct methods

- ✓ Direct Microscopic Count (Petroff-Hausser chamber).
- ✓ Standard plate count (viable count).
- ✓ Most probable number (MPN).
- ✓ Breed count method.
- ✓ Direct epifluorescent filter technique.
- ✓ Membrane filter technique.
- ✓ Flow cytometry.

Indirect methods


- ✓ Spectrophotometric method.
- ✓ Dry weight.
- ✓ ATP bioluminescence.
- ✓ Electrical impedance.

Standard plate count (viable count)

- A **viable** cell is defined as a cell which is able to divide and form a population (or colony).
- The standard plate count method consists of diluting a sample with sterile saline or phosphate buffer diluent until the bacteria are dilute enough to count accurately.
- **Serial dilution:** a set of dilutions in a mathematical sequence, as to obtain a culture plate with a countable number of separate colonies.
- The final plates in the series should have between 30 and 300 colonies.

Procedure

- Prepare eight test tubes each must contain 9ml of sterile normal saline.
- Transfer one ml of the original sample to the first tube by sterile pipette to form dilution number (1:10) or 10^{-1} (this is tenfold dilution).
- Transfer 1ml from 10^{-1} tube to second tube to form 10^{-2} dilution number, repeat the same step to remaining tubes to form dilutions 10^{-3} , 10^{-4} ,.....
- Take 0.1 ml from each dilution (shaking by vortex) and spread on the surface of several **Nutrient agar**, invert the plates & incubate at 37C° for 24 or 48 hours.



➤ Fewer than **30** colonies are not acceptable for statistical reasons, and more than **300** colonies on a plate are likely to produce colonies too close to each other to be distinguished as distinct **colony-forming units (CFUs)**.

➤ **Choosing plates that have colony between 25 and 250 colonies.**

➤ The number of bacteria in 1 ml of the original sample can be calculated as follows:

➤ Total No. of Bacteria (CFU)/ml = Average No. of colonies x 1/dilution or (Reciprocal of Dilution).

➤ **Suppose that you count 150 colonies in 10^{-5} dilution.**

$$\text{CFU/ml} = 150 \times \frac{1}{\text{dilution } (10^{-5})}$$

$$= 150 \times 10^5.$$

So CFU/1ml = 150×10^5 .

Colony counter



Some tips



- ✓ Some bacteria tend to stick together (two or more bacteria will give rise to one colony).
- ✓ If a sample has many different kinds of bacteria, it is not possible to have a medium or conditions that support various necessary growth conditions (blood agar & chocolate agar).
- ✓ Only viable organisms are counted.
- ✓ Only plates with colony numbers between 30-300 are useable.
- ✓ Use aseptic technique.
- ✓ Dilutions to be done accurately.

Direct Microscopic Count (Petroff-Hausser Chamber)



- In this method, a suspension of bacteria about 10 μl placed on a slide that has been ruled into squares and is designed to hold a specific volume of liquid.
- Sample could be diluted 1/10.
- By counting the bacteria that appear on the grid areas, the number of organisms in the sample can be calculated. It is a much faster test than the plate count.
- The drawback of this method is viable and nonviable organisms appear the same under a microscope.



- Chamber composed 9 big squares.
- Central big square used for bacterial counting which composed of 25 small squares each with 16 smaller ones.
- One big square is 1mm length, 1mm width.
- Space between the chamber to cover slide is 0.1mm.



Petroff-Hausser counting chamber



Cells calculation per 1 ml


- Suppose we count **150** cells in **5** small squares.
- Area of 1 large square = length \times width = $1\text{mm} \times 1\text{mm} = 1\text{mm}^2$
- 1 big square composed of 25 small squares, so area of 1 small = $1/25 = 0.04\text{mm}^2$.
- Volume of one small square = $0.04\text{mm}^2 \times 0.1\text{mm} = 0.004\text{mm}^3$.
- No of cells/ $1\text{mm}^3 = \text{total counted cells} / \text{no. of squares} \times \text{volume of 1 small square}$.
- No. of cells/ $1\text{mm}^3 = 150/5 \times 0.004 = 150/0.02 = 7500$ cells.
- To find the cells in 1cm^3 $\longrightarrow 1\text{cm}^3 = 1000 \text{ mm}^3$ so $7500 \times 1000 = 75 \times 10^5$ in 1cm^3
- So $1\text{cm}^3 = 1\text{ml}$
- No. of cells per $1\text{ml} = 75 \times 10^5/\text{ml}$.

Spectrophotometric (turbidity) method


- Is an instrument which measures an amount of light (photon) that a sample absorbs when it passes through the sample. To know the concentration of any substance.
- ✓ **Visible spectrophotometer: (400-700 nm).**
- ✓ **UV-visible spectrophotometer: (185-700 nm).**
- ✓ **IR spectrophotometer: (700 - 15000 nm).**
- Used to measure the **turbidity** or **optical density (O.D.)** of bacteria in a broth.
- This must be done with each bacteria because they are **different sizes**.
- sensitivity is limited to about 10^7 cells per ml for most bacteria.
- For instance, an optical density reading of **0.2 OD** for a broth culture of *E. coli* is nearly equal to **1 X 10^8 cells/ml**.
- For *S. aureus* **0.09 OD** is nearly equal to the same number.

Spectrometer

Photometer




Spectrophotometric (turbidity) method



Procedure

- ✓ Obtain an overnight fresh culture.
- ✓ Add 3 ml of an overnight culture to the first polystyrene cuvette.
- ✓ Add 3 ml of fresh media (without bacteria) to a second cuvette (used to zero the spectrophotometer).
- ✓ Turn it on, choose the wavelength of 600 nm button and enter.
- ✓ Place the cuvette with the fresh media into the sample chamber, read and zero it.
- ✓ Place your first cuvette with bacterial culture into the sampling chamber and read.
- ✓ Record the reading and use standard curve to know bacterial number .

Example, if you prepared a sample and got a spectrophotometer reading of 0.2 OD, you would find 0.2 on the Y-axis and read the bacterial count on the X-axis, in this case is equal to 1×10^8 cells/ml.



References

- ✓ Cappuccino, J.G. (2018). Microbiology A laboratory manual. 11th ed. Pearson Education, UK
- ✓ Madigan. M. T, Martinko. J. M, Stahl. D. A and Clark. D. P (2012) Brock Biology of Microorganisms. 13th ed. Pearson education, San Francisco.
- ✓ Schaechter. M (2010). The Desk Encyclopedia of Microbiology. Elsevier Ltd. California (USA).
- ✓ Benson (2001). Microbiological applications laboratory manual in General microbiology. 8th ed. The McGraw-Hill.
- ✓ <https://microbeonline.com/techniques-of-isolation-and-enumeration-of-bacteria/>
- ✓ <https://opentextbc.ca/microbiologyopenstax/chapter/how-microbes-grow/>
- ✓ <https://www.ruf.rice.edu/~bioslabs/methods/microscopy/cellcounting.html>
- ✓ https://wiki.hackuarium.ch/images/0/09/Bacterial_counts.pdf
- ✓ [https://chem.libretexts.org/Bookshelves/Physical_and_Theoretical_Chemistry_Textbook_Maps/Supplemental_Modules_\(Physical_and_Theoretical_Chemistry\)/Kinetics/02%3A_Reaction_Rates/2.01%3A_Experimental_Determination_of_Kinetics/2.1.05%3A_Spectrophotometry](https://chem.libretexts.org/Bookshelves/Physical_and_Theoretical_Chemistry_Textbook_Maps/Supplemental_Modules_(Physical_and_Theoretical_Chemistry)/Kinetics/02%3A_Reaction_Rates/2.01%3A_Experimental_Determination_of_Kinetics/2.1.05%3A_Spectrophotometry)
- ✓ <https://www.cwcboe.org/cms/lib/NJ01001185/Centricity/Domain/143/Molecular%20Biology/Intro%20to%20spectrophotometry/Introduction%20to%20Spectrophotometry.pdf>
- ✓ https://mathbench.umd.edu/modules-au/microbio_counting-methods/page04.htm#:~:text=A%20spectrophotometer%20doesn't%20exactly,get%20the%20total%20population%20size.