

Medical Gases

Prepared by: Dr. Azad J. Ali

Senior anesthesiologist

Gas supply: takes the form of either:

1-Cylinders




2-Piped gas system



Piped gas system

Bank of cylinders




Liquid Oxygen supply system



components of Medical gas cylinders

- 1. Removable metal safety cap.*
- 2. Cylinder valve.*
- 3. Safety pressure relief valve.*
- 4. Steel neck ring.*
- 5. Shoulders.*
- 6. Cylinder body.*



Removable metal cap



The body of the Cylinders

- ***Is made of seamless molybdenum steel.***
- ***It should be Rechecked every 5 years for testing tensile strength and the date of the check-up is marked with a label over the shoulder.***

The shoulder

Carries a label includes the following details :

- 1. Name ,chemical symbol , license , number and proportion of mixture.*
- 2. Batch number.*
- 3. Hazard warning and safety instructions.*
- 4. Cylinder size code.*
- 5. Cylinder content (liter).*
- 6. Maximum pressure (bars)*
- 7. Filling date , shelf life and expiry date.*
- 8. Directions for use.*
- 9. Storage and handling precaution.*

Cylinder label



The neck

- ***Ends in a tapered screw thread into which the valve is fitted .***
- ***The thread is sealed with a material that melt when the cylinder is exposed to intense heat, allowing the gas to escape and reducing the risk of explosion.***

The neck



Cylinder valve

- 1. Seals the cylinder content.***
- 2. The chemical formula of the particular gas is engraved on it.***
- 3. Have different types)the bullnose ,the hand wheel and the star) used in different circumstances.***



Pin-index
valve

Side spindle
pin-index
valve

Bullnose
valve



To Pin-index valve

Star valve

Valve types

Types of Cylinder valve

Components of gas cylinder valves

- 1. An On/off spindle to open and close the valve.***
- 2. Exit port for supplying gas to the medical apparatus (e.g. anesthetic machine).***
- 3. A safety relief valve to discharge excess contents to atmosphere if the cylinder is over pressurized.***
- 4. A non-interchangeable safety system (e.g. pin index system) is used on small size cylinders (E, F or even G) to prevent wrong gas connection.***

Pin index system

Pin index valves

Yoke on the anesthetic machine



Oxygen



Nitrous
oxide



Air



Entonox



The yoke



PROBLEMS IN PRACTICE AND SAFETY FEATURES OF CYLINDER VALVES

- 1. The exit port of cylinder valves may be occluded by dust , grease or dirt ,thus the plastic wrapping of the valve should me removed just before use or the valve must be slightly opened and closed before connecting to the anesthetic machine.***
- 2. The valve should be opened slowly to avoid rapid rise in the pressure and thus temperature of the gas in the machines pipelines, then fully opened (2 full revolutions).***
- 3. Avoid over tightening of the valve.***
- 4. Check Bodok seal for damage and have spare one.***

BODOK SEAL /WASHER FOR PIN INDEX REGULATOR



Cylinder sizes

Cylinders are manufactured in different sizes (from A to J).

Sizes A & H are not used for medical gases.

Size (E) cylinders are usually attached to the anesthetic machine.

- ***Size (J) cylinders are commonly used for cylinder manifolds.***
- ***Size (E) oxygen cylinders contain 680 L, whereas***
- ***Size (E) nitrous oxide release 1800 L.***
- ***Size (C) is the smallest size and can hold 1.2 L water.***
- ***Size (E) can hold 4.7 L ,while***
- ***Size (J) can hold 47.2L***



**Dimension
Hx/ Full wt.**

0.2m
2.3kg

0.5m
3.9kg

0.8m
6.5kg

1.3m
19kg

1.5m
78kg

**O₂ capacity at
177 bar**

200L

400L

600L

1400L

6000L

Water capacity

1.5L


2.9L

4.7L

23.6L

47.1L

Medical gas cylinders color codes



Medical gas cylinders new color codes

Example :





















- 1. Oxygen cylinders have White body and shoulder .***
- 2. Nitrous oxide cylinders have White body and Blue shoulder .***
- 3. CO2 cylinders have White body and Grey shoulder .***
- 4. Entonox (50% oxygen and 50% nitrous oxide) cylinders have White body and Blue/White quarters shoulder.***
- 5. Medical air cylinders have White body and White/Black quarters shoulder.***

British system

Color coding

GAS	SHOULDER	BODY
OXYGEN	WHITE	BLACK
NITROUS OXIDE	BLUE	BLUE
CYCLOPROPANE	ORANGE	ORANGE
CARBON DI OXIDE	GREY	GREY
AIR	WHITE	GREY
NITROGEN	BLACK	BLACK
ENTONOX	WHITE	BLUE

Color Coding of Gas Cylinders

Gas	USA	ISO*
Oxygen (O ₂)	 GREEN	 WHITE
Air	 YELLOW/SILVER	 WHITE/BLACK
Carbon dioxide (CO ₂)	 GRAY	 GRAY
CO ₂ and O ₂	 GRAY/GREEN	 GRAY/WHITE
Helium (He)	 BROWN	 BROWN
He and O ₂	 BROWN/GREEN	 BROWN/WHITE
Nitrous Oxide (N ₂ O)	 BLUE	 BLUE
Cyclopropane (C ₃ H ₆)	 ORANGE	 ORANGE
Ethylene (C ₂ H ₄)	 RED	 VIOLET
Nitrogen	 BLACK	 BLACK

Light weight cylinders

These are made from aluminum alloy with a fiberglass covering in epoxy resin matrix.

These can be used to provide O₂ at home, during transport or in MRI scanner.



Basic physics

- **Critical temperature**: of a substance is the temperature above which the substance can not be liquefied whatever pressure is applied. Critical temp. of Nitrous oxide and Oxygen are:
(36.5°C and -118°C) respectively.
- **Gas**: exists in the gaseous state at room temperature. Its impossible to change it into liquid at room temperature whatever pressure is applied since room temperature is above its critical temperature.
- **Vapor**: is the gaseous state of a substance below its critical temperature. At room temp. and atmospheric pressure the substance is liquid.
- **Filling ratio**: is weight of fluid in a cylinder divided by the weight of water required to fill the cylinder.

Cylinders content & delivery

- ***O₂ is stored as gas at a pressure of 13700 Kpa.***
- ***N₂O is stored a liquid with its vapor on top at a pressure of 4400 Kpa (the cylinder should be partially filled as liquid is less compressible.***
- ***Partial filling protects against risk of increase in pressure with increase ambient temp.***
- ***A full O₂ cylinder can deliver 130 times its capacity at atmospheric pressure.***
- ***Typical E size full O₂ cylinder delivering 4L/min will last for 2h &50 min , but at a rate of 15l/min can last only 45 min.***

At a constant temperature

- ***Gas containing cylinders shows a linear reduction in the pressure as it empties.***
- ***Liquid containing cylinders initially the pressure decrease constantly until the liquid part ends up, later on the pressure decreases.***
- ***The temperature in such cylinders decreases as a result of loss of latent heat of vaporization leading to ice formation at the outside of the cylinder.***

Problems in practice and safety features of medical gas cylinders

- ***Gases and vapors should be free of water vapor when stored in cylinders as it freezes at the exit ports on opening during cold weather.***
- ***color coded cylinders and pin index system prevents wrong gas connection to the anesthetic machine hence improving safety.***
- ***Cylinders must be checked regularly for sufficient content .***
- ***Cylinders must be stored in a purpose-built (dry ,well ventilated and fireproof)room, away from flammable materials as oil or grease ,they should not be exposed to continuous dampness , corrosive chemicals or fumes.***
- ***Full cylinders should be stored separately from empty ones to avoid mishaps.***
- ***Over pressurized cylinders are hazardous and must be reported to the manufacturer.***

Piped gas supply (piped medical gas and vacuum- PMGV)

- ***PMGV is a system where gases are delivered from central supply points to different sites in the hospitals at a pressure of about 400 Kpa through special outlet ports.***
- ***Oxygen ,nitrous oxide , Entonox, compressed air and vacuum are commonly supplied through the pipeline system.***

Piped gas supply

(piped medical gas and vacuum-PMGV)



Components

- 1. Central supply points (cylinder banks or liquid oxygen storage tanks).***
- 2. Network of pipelines made of special high-quality copper alloy(prevents degradation of the gases ,and has bacteriostatic properties. The fitting used are made of Brass and are brazed rather than soldering.***
- 3. The diameters of the pipelines differ according to the demand (42 mm for those leaving the manifold and 15 mm after repeated branching.***

Components

Bank of cylinders/manifold



Liquid oxygen



***pipelines made of special high-quality
copper alloy***



***42 mm for those leaving the manifold and 15
mm after repeated branching***

Components

4-Outlets are identified by color coding ,gas name and by shape (quick connect / disconnect probes) with an indexing collar specific for each gas (or gas mixture).

5-Outlet are installed as flush-fitting ,surface fitting units, on booms or pendants .

6-Flexible color-coded hoses connect to the outlets to the anesthetic machine. The machine end is fixed is thread specific non-interchangeable screw thread.

7-Area valve service units behind break glass covers at strategic points throughout the pipeline network. They can isolate the supply to an area in case of fire or an emergency.

(Quick connect / disconnect probes)





20220927_105033.mp4

(Quick connect / disconnect probes)



20220924_105222.mp4




20220924_112040.mp4



20220924_111856.mp4

Outlets are identified by color coding, gas name and by shape

outlets



probes



PMGV outlets

Booms




pendants



PMGV outlets

Booms



pendants






The machine end is fixed



The machine end is fixed



Area valve service units behind break glass covers



Area Valve Service Unit

Medical gas alarm panel

Main & area



Problems in practice and safety features

- 1. A reserve bank of cylinders is available if the primary supply fails.***
- 2. Single hose test for detection of cross connection.***
- 3. Tug test to detect misconnection***
- 4. Regulations for PMGV installation , repair and modification are enforced.***

Tug test

- Connect O₂ pipeline to oxygen wall outlet using the Schrader quick coupler system
- Correct coupling will not allow detachment of pipeline from the Schrader coupler when a tug is given to the pipeline.
- Similar test can be performed with N₂O pipeline with N₂O wall outlet.

Problems in practice and safety features

5- Anesthetists are responsible for gases supplied from the terminal outlet through the anesthetic machine. Pharmacy and engineering departments share the responsibility for gas pipelines 'behind the wall'.

6-There is a risk of fire from worn or damaged flexible high pressure hoses that transfer gases from the cylinders or outlets to the ventilators or anesthetic machine, regular inspection and replacement , every 2-5 years, of all gases hoses.

Sources of gas supply

Sources in central pipeline system can be:

- 1. Cylinder manifold(s).including other gases***
- 2. Liquid oxygen storage tank.***
- 3. Oxygen concentrator.***

Cylinder manifold components

- 1. Two equal groups of J-size cylinders (primary and secondary), alternate in supplying the pipeline. The number depends on the demand.***
- 2. Cylinders in each group are connected to a common pipe through a non-return valve, the two common pipes are connected through a pressure regulator to the pipeline.***
- 3. Nitrous oxide is only available in cylinders thus its manifold is larger. O2 manifold acts as a back up for liquid O2 supply.***

Mechanism of action

- 1. All the cylinders in either group are opened ,this allows them to empty simultaneously.***
- 2. When the primary group is nearly empty, the supply is automatically changed to the secondary group through a pressure sensitive device.***
- 3. The change over activates an electrical signaling system to alert the staff that cylinders change is needed.***

Problems in practice and safety features

- 1. Manifold should be housed in well-ventilated and fire proof building away from the hospital building.***
- 2. Manifold room should not be used as a general cylinder store.***
- 3. All empty cylinder should be removed immediately from the manifold room.***

Liquid Oxygen

A vacuum insulated evaporator (VIE) or giant thermos flask is the most economical way to store and supply O₂.



Components:

- 1. A thermally insulated double walled steel tank with a layer of perlite in a vacuum is used as the insulator.*
- 2. A pressure regulator maintains the pipeline pressure of about 400 kpa.*
- 3. A safety valve that opens at 700 kpa relieves build-up pressure within the vessel which results from under demand for O₂.*
- 4. A control valve opens when there is extra demand on O₂. this allows liquid O₂ to evaporate by passing through super heater made of insulated coils of copper tubing.*

Perlite



Liquid Oxygen



Mechanism of action

- 1. Liquid O₂ is stored (up to 1500 L) at (-150 to -170)C lower than the critical temperature and at a pressure of 10.5 bars.***
- 2. The temp. of the vessel is maintained by the high vacuum shell. Evaporation of liquid O₂ needs heat (latent heat of vaporization) which is taken from the liquid O₂ helping to maintain its low temp.***
- 3. The contents of the vessel can be calculated through a special device. When required fresh liquid O₂ can be pumped from a tanker to the vessel.***
- 4. The cold O₂ gas is warmed in a coiled copper tubing, once outside the vessel, the rise in the temp. causes an increase in the pressure.***
- 5. At a temp. of 15C and atmospheric pressure ,liquid O₂ can give 842 times its volume as gas.***

Problems in practice and safety features

- 1. Reserve banks of cylinders are kept in case of supply failure.***
- 2. A liquid O₂ storage vessel should be away from the hospitals main building.***
- 3. Spillage of cryogenic liquid can cause cold burns or frostbite and hypothermia.***

Oxygen concentrators

- ***Extract O₂ from air by differential adsorption.***
- ***Either small that supply O₂ for a single patient or anesthetic machine or ;***
- ***Large to supply O₂ for a medical gas pipeline system.***

Components:

many columns of hydrated Aluminum silicates of the alkaline earth metals in a powder or granular form.


Small O₂ concentrator single patient use



Large O₂ concentrator



Large O₂ concentrator



Mechanism of action

- 1. Ambient air is pressurized by a compressor at 137 kpa .***
- 2. Air is exposed to the filter that is composed of the aluminum columns ,it retains nitrogen and the other unwanted components that are released to the atmosphere .***
- 3. The maximum O₂ concentration is 95%.***
- 4. The life of the filter is at least 20 000 hours or 10 years.***
- 5. Routine maintenance is to change the filter.***

Entonox

- ***Is a compressed gas mixture of 50% O₂ and 50% N₂O.***
- ***Its commonly used in the casualty and labor room to provide analgesia.***
- ***A 2-stage pressure demand regulator is attached to it.***
- ***The gas flow =s only when the patient inspires through a mask or mouth piece.***
- ***Entonox is compressed at 13700 kpa.***
- ***The cylinders should be stored at 10`C for 24 hours before use.***

Entonox



Entonox

Mouth piece



Face mask



Problems in practice and safety of Entonox

- ***At temperatures below -5.5 C liquefaction and separation of the two components occur, forming a vapor part rich in O₂ and a liquid part with poor O₂ content, thus rewarming and shaking reverses separation and liquefaction.***
- ***The horizontal positioning increases the surface area for diffusion.***
- ***If the contents are mixed by repeated inversion it can be used earlier than 24 hours.***

Compressed air

- ***Medical air is supplied in hospitals for clinical use at a pressure of 400 kpa***
- ***or to drive tools like orthopedic saw or drill at a pressure of 700 kpa***
- ***The terminal outlets for each is different to avoid misconnection.***
- ***Medical air is oil free ,cleaned by special filters and separators.***
- ***Its supplied from cylinder manifolds, or more economic compressor plant with duty and back-up compressors.***

Compressed Medical air

